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**Cat:** G1366

**Size:** 3×20mL /3×50mL

**Storage:** 2-8℃, avoid light, valid for 6 months.

### Kit Components

Reagent	3×20mL	3×50mL	Storage
Reagent(A): Oxidant	20mL	50mL	2-8℃, avoid light
Reagent(B): Schiff Reagent	20mL	50mL	2-8℃, avoid light
Reagent(C): Sodium Sulfite Solution	20mL	50mL	RT, avoid light

### Introduction

Glycogen Staining is one of the routine staining methods in pathology. It is often used to display glycogen and other polysaccharides, neutral mucous substances and some acidic substances, cartilage, pituitary, fungi, pigments, amyloid substances, basement membrane and so on.

The dyeing principle of Glycogen Periodic Acid Schiff (PAS) Stain Kit (For Fungus) is that oxidant oxidizes polysaccharides on fungi wall and exposes aldehyde groups, then aldehyde groups combine with Schiff Reagent and show red color. This kit has stable performance and strong specificity, simple operation, only for 0.5 h, low concentration oxidant is more suitable for fungal staining, no hydrochloric acid ethanol differentiation solution differentiation steps.

### Self Provided Materials

10% Neutral formalin, Distilled water, Series of ethanol

### Protocol (for reference only)

1. Conventional fixation (usually 10% formalin), dehydration and embedding.
2. Drop the cell smear into distilled water (without embedding or dewaxing) or dewax the tissue section in distilled water.
3. Oxidize in Oxidant at room temperature for 5-8mins.
4. Rinse with tap water twice and soak with distilled water twice.
5. Add Schiff Reagent and cover the slice. Place it in dark place at room temperature for 10-20mins.
6. Drip with Sodium Sulfite Solution twice for 2mins each time.
7. Rinse with running water for 5mins.(according to the microscopic observation)
8. Dehydrate with series of ethanol, transparent by xylene and seal with resinene.

### Result

Fungus	Purplish Red
Cytoplasm	Red in different degrees

### Negative Control (Optional)

1. Drop 2mL G1284-Diastase Solution, 1%, pH5.3, Water Solvent and treat 37℃ for 60mins, then oxidize with other samples. The result should be negative.
2. (Alternative) Treat saliva tablets (filtered and used) for 30-60mins, then oxidize with other samples. The result should be negative.

### Note

1. The oxidation time should not be too long. The optimum temperature for oxidation is 18-22℃.
2. Avoid excessive exposure to sunlight and air when using Oxidant and Schiff Reagent. It is better to remove and restore to room temperature 30mins before use and use in dark places.
3. The action time of Oxidant and Schiff Reagent is very important, which depends on the thickness of section, the type of cell or tissue, etc.
4. For routine sections, Glycogen PAS Staining Solution is recommended because of its relatively high concentration of Oxidant and Mayer Hematoxylin Staining Solution.
5. For your safety and health, please wear experimental clothes and disposable gloves.