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Tartrate-Resistant Acid Phosphatase (TRAP) Stain Kit

Cat: G1492

Size: 4×10mL/4×20mL

Storage: -20°C, avoid light, valid for 3 months.

Kit Components

Reagent		4×10mL	4×20mL	Storage
Reagent(A): TRAP Fixative		50mL	100mL	2-8°C, avoid light
Reagent(B): TRAP Incubation Solution	B1: AS-BI Buffer	1mL	2×1mL	-20°C, avoid light
	B2: GBC Solution	0.1mL	0.2mL	-20°C, avoid light
	B3: TRAP Buffer	9mL	18mL	RT, avoid light
Mix reagent B1,B2,B3 in 10:1:90 ratio as TRAP Incubation Solution before use and not store for long.				
Reagent(C): Hematoxylin Solution		10mL	20mL	2-8°C, avoid light
Reagent(D): Methyl Green Solution		10mL	20mL	RT, avoid light
Reagents C and D are both re staining solutions, and one can be used during the staining process. Repeated staining is not recommended.				

Introduction

Acid phosphatase (ACP) is widely distributed in various tissues, mainly in the lysosome, so it is often used as a lysosome marker enzyme. The acid phosphatase outside of lysosome exists in endoplasmic reticulum and cytoplasm. Acidic phosphatase is different in all kinds of animals. The optimum pH for its activity is 4.5-5.5. Tartrate-resistant acid phosphatase (TRAP), which mainly exists in normal human alveolar macrophages and leukemia human spleen cells, is released into the cell follicles rather than into the blood. The majority of TRAP in blood comes from osteoclasts, so we can know the function of osteoclasts by measuring TRAP in blood.

The reaction principle of Tartrate-Resistant Acid Phosphatase (TRAP) Stain Kit is basically the same as that of acid phosphatase staining. On the basis of conventional staining, tartrate resistant acid phosphatase was introduced for enzyme subtype screening. Only tartrate resistant acid phosphatase active sites were stained on cell or tissue sections, while conventional acid phosphatase sites were not stained. This staining solution can be used for fresh blood and cell smears, as well as for frozen sections and paraffin sections.

Self Provided Materials

Distilled Water, Incubator, Slide, Microscope

Protocol(for reference only)

Sample Processing

1. Cell smear: Take fresh whole blood or bone marrow to prepare a smear, place the slide at 30 degrees on the slide, directly in front of the blood or cell droplet, move it back slightly to contact the blood or cell droplet, and make the latter scatter along the lower edge of the slide. Then, slide it smoothly forward along the slide plane until evenly spread, and let it dry naturally.
2. Adherent culture of cells and cell creep: remove the culture medium, wash with 1×PBS 3-4 times, and drain excess water. Remove the cell slides.
3. Frozen slice: Remove from the refrigerator and place at room temperature. Soak in distilled water and wash for 1-2 min, then drain excess water.
4. Paraffin section: Paraffin sectioning is dewaxed for 5-10 min and repeated once. Anhydrous ethanol for 5 min, 90% ethanol and 70% ethanol for 2 min each, soak in distilled water for 2 min, and drain excess water.

Staining Incubation

1. Fix with pre-cooled TRAP Fixative at 2-8 °C for 30s-3min, and in most cases, 30-60s is sufficient.
2. Wash with distilled water and drain excess water (not too dry).
3. (Slicing treatment) Drip TRAP incubation solution onto the slices and incubate them in a wet box at 37 °C for 45-60 min. Wash with distilled water.
4. (Adherent cells) Drip TRAP incubation solution to cover the cells, incubate at 37 °C for 45-60 min, and wash with distilled water.
5. Re staining: Stain with Hematoxylin Solution or Methyl Green Solution for 2-3 min. If Hematoxylin Solution is used for re staining, must return blue with tap water for 10 min or with 1×PBS for 3-5 min after

staining.

Sealing observation

1. (Slicing treatment) Wash slightly with distilled water, drain excess water, add water-based sealing agent dropwise, seal the film, and observe under a microscope.
2. (Adherent cells) Wash slightly with distilled water and observe under a microscope with water or 1×PBS.

Result

Positive	Purplish Red
Nucleus	Blue(Hematoxylin) or Green(Methyl Green)

Clinical Significance

1. Lung resident macrophages are usually TRAP positive, and their frozen or paraffin sections can be used as positive control sections.
2. Multiple hairy cells in peripheral blood of individuals with polycystic leukemia showed TRAP positivity.
3. In the bone samples fixed and decalcified by enzyme protection, osteoclasts and chondroclasts showed TRAP positivity.

Note

1. This staining kit is an enzyme active in situ staining kit commonly used for cell identification and specific cell localization in tissues after osteoclast induction culture. The reaction substrate and dye are prone to failure, and it is recommended to store the product properly according to the temperature described in the instruction manual as soon as received.
2. Reagent B2 is a reserve solution for azo dyes, and the presence of granular suspensions is a normal phenomenon that does not affect dyeing. If you are concerned, pre mixing B1 and B2 can completely dissolve the particulate matter, but it needs to be used within 2 hours.
3. Enzyme activity is easily affected by various factors, leading to attenuation or disappearance. It is recommended to take fresh samples and use pre cooled fixing solution to fix them in a refrigerator at 2-8 °C for no more than 24 hours. It is recommended to use low melting point wax (52-56 °C) during embedding to avoid thermal deactivation.
4. For your safety and health, please wear laboratory clothes and disposable gloves for operation.