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Weigert's Gram Stain Kit(Conn Modified)

Cat: G3220

Size: 4×50mL

Storage: RT, avoid light, valid for 6 months.

Kit components

Reagent	4×50mL	Storage
Reagent(A): Carmine Staining Solution	50mL	RT, avoid light
Reagent(B): Acid Alcohol Differentiation Solution	50mL	RT
Reagent(C): Ammonium Oxalate Crystal Violet Staining Solution	50mL	RT, avoid light
Reagent(D): Weigert Iodine Solution	50mL	RT, avoid light

Introduction

Gram staining was invented by Danish doctor Christian Gram in 1884. It is not only a differential staining method widely used in bacteriology, but also a double staining method. The bacteria can be identified as Gram-positive bacteria and Gram-negative bacteria by this staining method. The different color reaction of bacteria is due to the difference of cell wall permeability to ethanol and anti decolorization ability, which is mainly determined by the thickness and structure of peptidoglycan layer. In the same dyeing environment, using different isoelectric points of bacteria (the isoelectric points of gram positive bacteria are pH 2.0-3.0, and the isoelectric points of gram negative bacteria are pH 4.0-5.0), the negative charge of positive bacteria is more than that of negative bacteria, which is firmly combined with positively charged alkaline dyes such as crystal violet. After adding mordant (iodine) into the bacteria, It combines with dyes to form water-insoluble crystal violet iodine protein complex, and combines with ribonucleic acid magnesium salt in positive bacteria, so that colored bacteria are not easy to decolorize; The differentiation solution(aniline, acetone, etc.) is not easy to penetrate the cell wall of the positive bacteria, so the positive bacteria are not easy to fade, but the differentiation agent is easy to enter the negative bacteria to dissolve the dye and iodine complex to decolorize the negative bacteria.

Weigert's Gram Stain Kit(Conn Modified) is improved from the classical gram dyeing formula, using carmine as the lining agent. This method is also evolved from the ammonium oxalate crystal violet method. The clinical specimens can be smeared directly, the background is clean, the intracellular phagocytes are clear and easy to identify, and the bacterial staining characteristics are typical. Gram positive bacteria and gram negative bacteria can be distinguished, especially suitable for identifying bacterial and non bacterial blue micro granular substances (such as calcium salt).

Self Provided Materials

10% formalin fixative, optical microscope.

Protocols(for reference only)

1. Fix the tissue in 10% formalin fixative, dehydrate and embed by routine method.
2. Cut slice in thickness of 4-5 μ m. Conventionally dewax to water.
3. Drip Carmine Staining Solution onto the section and stain for 5min, then remove the solution.
4. Treat with Acid Alcohol Differentiation Solution for 2-5s.
5. Rinse with distilled water for 3-5s.
6. Drip dyeing with Ammonium Oxalate Crystal Violet Staining Solution for 5 min, and remove the dyeing solution.Slightly suck dry the remaining solution around the slice with filter paper.
7. Drip Weigert Iodine Solution onto the slice for 1min, remove the iodine solution and repeatedly suck dry the water on the slice with filter paper.
8. Decolorize with aniline xylene solution (1:2), shake gently from time to time until there is no purple prolapse on the slice, wash with xylene immediately and observe under the microscope.
9. Drop in new xylene, wash repeatedly for many times, wash away aniline thoroughly, and seal with resinene.

Result

Gram positive bacteria and cellulose	Blue Purple
Gram negative bacteria	Red

Nucleus

Red

Note

1. After staining with Carmine Staining Solution, without washing, directly differentiate with Acid Alcohol Differentiation Solution to avoid fading.
2. When differentiating with aniline xylene solution, gently shake the slices to make the differentiation uniform. If the differentiation is slow, can replace new aniline xylene. Until there is no purple solution out of the slice, immediately pour aniline xylene and drip new xylene for washing.
3. After washing with xylene, shall observe under the microscope. In case of insufficient differentiation, aniline and xylene can be added to continue differentiation until Gram-positive bacteria are clear, but be careful not to over differentiate.
4. Finally, wash the slices repeatedly with xylene to completely remove the aniline. If a small amount of aniline remains in the slice, it will fade easily in the future.
5. Gram positive bacteria must have an undamaged bacterial wall. If the bacterial wall is damaged, the staining is negative.
6. For your safety and health, please wear experimental clothes and disposable gloves.