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Alizarin Red S Staining Quantitative Detection Kit (CPC Method)

Cat: G3283

Size: 2×20mL/2×100mL

Storage: 2-8°C, avoid Light, valid for 6 months.

Kit Components

Reagent	2×20mL	2×100mL	Storage
Reagent(A): Stain Solution	20mL	100mL	2-8°C, avoid Light
Reagent(B): Standard Solution	1mL	5mL	2-8°C, avoid Light
Reagent(C): Elution Solution	20mL	100mL	RT

Introduction

Alizarin Red S is an anthraquinone derivative, which is a sodium salt of Alizarin Sulfonate. It can chelate with calcium salts in calcium carbonate or calcium phosphate to form an orange red complex, suitable for tissue sectioning and calcium salt deposition staining of cultured cells. For the convenience of analyzing calcium salt content and concentration, Alizarin Red S staining can be completely eluted with CPC reagent. The eluted solution is placed in a 96 well plate, and the calcium salt content is semi quantitatively analyzed by measuring the OD value at 560nm.

This reagent kit consists of Stain Solution, Standard Solution, and Elution Solution, and it is suitable for Alizarin Red S staining, extraction and determination experiments of calcium salts in cell and tissue samples.

Protocol(for reference only)

Reagent(C): Elution Solution is prone to crystallization in low temperature environments. It is recommended to store it separately at room temperature. If precipitation is observed, it can be dissolved in a 37 °C water bath or a 60 °C oven after heating.

1. For tissue sections: paraffin section thickness 4 μm. Conventional dewaxing to water. For cultured cells, they are fixed at room temperature for 15-30 minutes using 4% paraformaldehyde or 10% neutral formalin fixative. Rinse with distilled water twice for 3min each.
2. Add Reagent(A): Stain Solution and stain for 3-5 minutes. Then quickly wash off excess staining solution.
3. For tissue sections, add 400ul Reagent(C): Elution Solution onto the section and incubate at room temperature for 15-30 minutes; for cultured cells, six well plate add 1ml Reagent(C): Elution Solution and incubate at room temperature for 15-30 minutes, until the Alizarin Red S staining is completely washed off.
4. Collect the eluted Alizarin Red S solution using a centrifuge tube and shake well. Place in a 96 well plate, add 100uL to each well, with a total of three repeated wells, and measure the OD value at a wavelength of 560nm on an enzyme-linked immunosorbent assay (ELISA) reader.

Standard curve preparation

Dilute Reagent(B): Standard Solution by 20 times with Reagent(C): Elution Solution to prepare standard solution of 0.5g/L), and then perform multiple dilution in sequence to obtain standard solutions of 0.25g/L, 0.125g/L, 0.0625g/L, 0.03125g/L and 0.015625g/L. Use Reagent(C): Elution Solution as the blank control. Place in a 96 well plate, add 100uL to each well, measure the OD value at a wavelength of 560nm on the enzyme-linked immunosorbent assay (ELISA), and organize the data to develop a standard curve.

Result

Substitute the OD value obtained from the test sample into the standard curve linear equation to obtain the corresponding Alizarin Red S concentration of the test sample, thereby semi quantitatively determining the calcium salt content.

Note

1. The determination experiment of the same batch of samples should be conducted on the same day as much as possible. For samples with an interval of more than three days, it is best to redefine the standard curve.
2. For your safety and health, please wear laboratory clothes and disposable gloves when operating.

Related products

G1450 Alizarin Red S Solution, 0.2%, pH8.3

G1452 Alizarin Red S Solution, 1%, pH4.2

G1453 Alizarin Red S Staining Extraction Solution (CPC, 10%)

G3280 Calcium Stain Kit (Modified Alizarin Red S Method)

G3281 Calcium Stain Kit (Alizarin Red S Method)