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规格：3×10mL/3×100mL

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Storage: RT, valid for 1 year.

Kit Components

| Reagent | 3×10mL | 3×100mL | Storage |
|------------------------------------|--------|---------|---------|
| Reagent(A): Blue marking solution | 10mL | 100mL | RT |
| Reagent(B): Green marking solution | 10mL | 100mL | RT |
| Reagent(C): Black marking solution | 10mL | 100mL | RT |

Introduction

The cut surface exposed after tissue removal is called the cutting edge, which is commonly used in tumor resection surgery. When the tissue is first cut off, its shape remains intact and the orientation of the cutting edge can be clearly distinguished. However, after fixation, gradual dehydration, and slicing, it becomes more difficult to distinguish the orientation of the cutting edge. However, when conducting macroscopic examination of the tissue removed by surgery, it is usually necessary to permanently determine the direction of the specimen for further macroscopic examination and microscopic identification of the resection edge.

This kit is a combination of three labeled dyes that can be clearly contrasted with tissue under a microscope, namely blue, green, and black. Moreover, after being dehydrated and embedded, the tissue labeled with this reagent still has sufficient color retention and high staining fastness. This product can choose to directly apply markings on fresh tissue, or mark after fixing the tissue. Choosing to mark after fixing the tissue can better ensure the stability of the organizational structure.

Protocols(for reference only)

For unfixed sample

1. Wash the fresh tissue obtained with physiological saline to remove residual blood on the cut surface.
2. Use filter paper or tissue to absorb residual moisture on the surface, and if necessary, wipe the cut surface with acetone or 2% glacial acetic acid to improve the adhesion of the marking solution.
3. Use a cotton swab to dip an appropriate amount of the selected color marking solution onto the cutting edge, and let it stand for 3-5 min to fully adhere the dye into a film.
4. Liquid nitrogen or dry ice ethanol mixed slurry for rapid freezing to prepare frozen slices.

For fixed sample

1. After the tissue sample is fully fixed, remove it and use filter paper or tissue to absorb residual surface moisture. If necessary, acetone or 2% glacial acetic acid can be used to wipe the cut surface to improve the adhesion of the labeling solution.
2. Select the marking solution of the selected color to apply or dip it onto the tissue cutting edge, and let it stand for 3-5 min to fully adhere the dye into a film.
3. Place the labeled tissue back into the fixative and fix for at least 1 h.
4. You can choose to continue fixing the tissue or directly dehydrate and embed it.

Note

1. Acetone or 2% glacial acetic acid can enhance the adhesion of the labeling solution, but it can also cause slight dehydration and deformation of the edges. It can be used according to the sample condition.
2. The thickness of the marking pigment layer should be appropriate to ensure sufficient coloring, but it should not be too thick to affect the drying speed.
3. The labeled tissue should avoid direct contact with the labeled coating using forceps or other hard objects as much as possible.
4. For your safety and health, please wear laboratory clothes and disposable gloves to exercise.